

Anti-Nod1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y239

Item Number: CM0013271



Introduction

High-quality rabbit polyclonal antibody against human NOD1 (CARD4), a pattern recognition receptor for bacterial peptidoglycan. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. For innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	NOD1 Rabbit pAb
Remarks/Alias	CARD4; NLRC1; CLR7.1; NOD1
Species	Human
Gene ID	10392
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human NOD1 (NP_006083.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y239
Protein Weight	108kDa
Shipping	Ice bag