


Anti-Nppa Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01160

Item Number: CM0014630



Introduction

High-quality rabbit polyclonal antibody against human NPPA (Natriuretic peptides A), validated for WB, IHC-P, ELISA. Reacts with mouse and rat. Ideal for cardiovascular and renal research.

[Learn More](#)

Parameter	Value
Product Name	NPPA Rabbit pAb
Remarks/Alias	ANF; ANP; CDD; CDP; PND; ATFB6; ATRST2; CDD-ANF; NPPA
Species	Human
GeneID (human)	4878
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NPPA (NP_006163.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01160
Protein Weight	16 kDa
Shipping	Ice bag