

Anti-Nr0B1 Polyclonal Antibody For Wb, If/Icc, Elisa - P51843

Item Number: CM0024477



Introduction

Rabbit polyclonal antibody targeting human NR0B1 (DAX1), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for studying nuclear receptor corepression and HPA axis development.

[Learn More](#)

Parameter	Value
Product Name	NR0B1 Rabbit pAb
Remarks/Alias	AHC; AHX; DSS; GTD; HHG; AHCH; DAX1; DAX-1; NROB1; SRXY2; NR0B1
Species	Human
GeneID (Human)	190
GeneID	190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human NR0B1 (NP_000466.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51843
Protein Weight	52kDa
Shipping	Ice bag