


Anti-Nrp2 Polyclonal Antibody For Wb, If/Icc, Elisa - O60462

Item Number: CM0029180



Introduction

Rabbit polyclonal antibody against human Neuropilin-2 (NRP2/VEGF165R2). Suitable for Western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat targets. For vascular and neuronal research applications.

[Learn More](#)

Parameter	Value
Product Name	NRP2 Rabbit pAb
Remarks/Alias	NP2, NPN2, PRO2714, VEGF165R2, NRP2
Species	Human
Gene ID (Human)	8828
Gene ID	8828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-219 of human NRP2 (NP_957718.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60462
Protein Weight	105kDa
Shipping	Ice bag