


Anti-Nsd2/Whsc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O96028

Item Number: CM0030375



Introduction

Rabbit polyclonal antibody against WHSC1/NSD2, a histone methyltransferase. Suitable for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse. For epigenetic and chromatin biology studies.

[Learn More](#)

Parameter	Value
Product Name	WHSC1 Rabbit pAb
Remarks/Alias	WHS; TRX5; KMT3F; KMT3G; MMSET; RAUST; WHSC1; REIIBP
Species	Human
GeneID (Human)	7468
GeneID	7468
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human WHSC1 (NP_579877.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O96028
Protein Weight	152kDa
Shipping	Ice bag