


Anti-Nsdhl Monoclonal Antibody For Wb, Elisa - Q15738

Item Number: CM0003200



Introduction

Research-grade anti-NSDHL rabbit monoclonal antibody validated for WB and ELISA. Targets human NSDHL (Q15738), a sterol-4-alpha-carboxylate 3-dehydrogenase involved in cholesterol biosynthesis and EGFR trafficking.

[Learn More](#)

Parameter	Value
Product Name	NSDHL Rabbit mAb
Remarks/Alias	H105E3; XAP104; SDR31E1; NSDHL
Species	Human
GeneID	50814
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human NSDHL (NP_057006.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15738
Protein Weight	42kDa
Shipping	Ice bag