

Anti-Nsun2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q08J23

Item Number: CM0028055

Introduction

High-quality NSUN2 Rabbit pAb validated for WB, IHC-P, IF/ICC, and ELISA. Targets human NSUN2 (Q08J23) with cross-reactivity to mouse and rat. Ideal for RNA methylation and differentiation research.

[Learn More](#)

Parameter	Value
Product Name	NSUN2 Rabbit pAb
Remarks/Alias	MISU; MRT5; SAKI; TRM4; NSUN2
Species	Human
GeneID (Human)	54888
GeneID	54888
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 617-708 of human NSUN2 (NP_060225.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08J23
Protein Weight	86kDa
Shipping	Ice bag