


Anti-Nth1 Rabbit Monoclonal Antibody For Wb, Elisa - P78549

Item Number: CM0003171



Introduction

Recombinant rabbit monoclonal antibody against human NTH1 (Endonuclease III-like protein 1) for Western blot and ELISA applications. Validated in human samples, suitable for DNA repair and oxidative damage research.

[Learn More](#)

Parameter	Value
Product Name	NTH1 Rabbit mAb
Remarks/Alias	FAP3; NTH1; OCTS3; hNTH1
Species	Human
GeneID (Human)	4913
GeneID	4913
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-136 of human NTH1 (NP_002519.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78549
Protein Weight	34kDa
Shipping	Ice bag