

Anti-OCrl Rabbit Polyclonal Antibody For Wb, Elisa - Q01968

Item Number: CM0020931



Introduction

High-quality rabbit polyclonal antibody targeting human OCRL (Q01968). Validated for WB and ELISA. Cross-reacts with mouse and rat. Detects the inositol polyphosphate 5-phosphatase involved in endosomal trafficking and ciliary function. Suitable for Lowe syndrome research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] OCRL Rabbit pAb
Remarks/Alias	LOCRL; DENT2; NPHL2; OCRL1; Dent-2; INPP5F; OCRL-1; RL
Species	Human
Gene ID (Human)	4952
Gene ID	4952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human OCRL (NP_001578.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01968
Protein Weight	104kDa
Shipping	Ice bag