


Anti-Ogg1 Rabbit Polyclonal Antibody For Wb And Elisa - O15527

Item Number: CM0020808



Introduction

Rabbit polyclonal antibody against human OGG1 (N-glycosylase/DNA lyase) for DNA repair research. Validated for WB and ELISA, reacts with human and mouse OGG1. Ideal for oxidative damage studies.

[Learn More](#)

Parameter	Value
Product Name	OGG1 Rabbit pAb
Remarks/Alias	HMMH; MUTM; OGH1; HOGG1; OGG1
Species	Human
GeneID (human)	4968
GeneID	4968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-345 of human OGG1 (NP_002533.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15527
Protein Weight	39kDa
Shipping	Ice bag