

Anti-Olig2 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q13516

Item Number: CM0010971

Introduction

Rabbit monoclonal antibody against human OLIG2, validated for Western blot, immunofluorescence, IHC, and ELISA. Suitable for studying oligodendrocyte and motor neuron specification, gliogenesis, and gliomagenesis. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	Olig2 Rabbit mAb
Remarks/Alias	BHLHB1; OLIGO2; RACK17; PRKCBP2; bHLHe19
Species	Human
GeneID (Human)	10215
GeneID	10215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 168-323 of human OLIG2 (NP_005797.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13516
Protein Weight	32kDa
Shipping	Ice bag