

Anti-OsdlT Polyclonal Antibody For Wb And Elisa - Q9Lwu9

Item Number: CM0020577

Introduction

This rabbit polyclonal antibody specifically recognizes the OsDLT protein (SWISS: Q9LWU9) in rice. Validated for WB and ELISA applications, it is an essential tool for studying brassinosteroid signaling, transcription regulation, and plant organ size control.

[Learn More](#)

Parameter	Value
Product Name	OsDLT Rabbit pAb
Remarks/Alias	DLT; OsGRAS-32; D62; GS6; SMOS2; OsDLT
Species	Oryza sativa
Gene ID	4339978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-214 of Oryza sativa OsDLT (Q9LWU9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q9LWU9
Protein Weight	66kDa
Shipping	Ice bag