

Anti-Oseil1 Rabbit Polyclonal Antibody For Wb, Elisa - Q10M41

Item Number: CM0024450

Introduction

OsEIL1 Rabbit Polyclonal Antibody validated for WB and ELISA in rice (*Oryza sativa*). Ideal for studying ethylene signaling, stress response, and disease resistance pathways. Target MW 71kDa.

[Learn More](#)

Parameter	Value
Product Name	OsEIL1 Rabbit pAb
Remarks/Alias	EIL1; MHZ6; EIL1A; OsEIL1; OsEIL1a; OsJ_10672
Species	<i>Oryza sativa</i>
GeneID (human)	EIL1
GeneID	107276358
Immunogen	Recombinant protein LSRGKDPAGGVVGDPSC-C
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q10M41
Protein Weight	71kDa
Shipping	Ice bag