

Anti-Osrd6/Shl2 Rabbit Polyclonal Antibody For Wb & Elisa - Q8Lhh9

Item Number: CM0024182

Introduction

Rabbit polyclonal antibody targeting rice OsRDR6 (SHL2, SHOOTLESS2) for Western blot (WB) and ELISA. This RNA-dependent RNA polymerase is critical for shoot meristem development and the ta-siRNA pathway, making it suitable for plant developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	OsRDR6 Rabbit pAb
Remarks/Alias	SHL2; SHOOTLESS2; OsRDR6
Species	Oryza sativa
Gene ID	4323938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 963-1218 of Oryza sativa OsRDR6 (Q8LHH9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q8LHH9
Protein Weight	137kDa
Shipping	Ice bag