


Anti-Otc Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P00480

Item Number: CM0005091



Introduction

Recombinant rabbit monoclonal OTC antibody validated for WB, IF-P, and ELISA with cross-reactivity to human, mouse, and rat. Essential for urea cycle research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	OTC Rabbit mAb
Remarks/Alias	OCTD; OTC1; OTCD; OTCas; OTC
Species	Human
Gene ID	5009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 215-354 of human OTC (NP_000522.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00480
Protein Weight	40kDa
Shipping	Ice bag