

Anti-Otx2 Rabbit Polyclonal Antibody For Wb And Elisa - P32243

Item Number: CM0028067

Introduction

OTX2 Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reactive with mouse. Targets human homeobox protein OTX2 (32 kDa), ideal for developmental neuroscience and rare disease research.

[Learn More](#)

Parameter	Value
Product Name	OTX2 Rabbit pAb
Remarks/Alias	CPHD6; MCOP55; OTX2
Species	Human
Gene ID (Human)	5015
Gene ID	5015
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 140-280 of human OTX2 (NP_068374.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P32243
Protein Weight	32kDa
Shipping	Ice bag