


Anti-Oxa1L Rabbit Polyclonal Antibody For Wb And Elisa - Q15070

Item Number: CM0027217



Introduction

Anti-OXA1L Rabbit Polyclonal Antibody (pAb) specifically detects human OXA1L and cross-reacts with mouse and rat. Validated for Western blot (WB) and ELISA, ideal for mitochondrial inner membrane protein insertase research.

[Learn More](#)

Parameter	Value
Product Name	OXA1L Rabbit pAb
Remarks/Alias	OXA1; OXA1L1; OXA1L
Species	Human
Gene ID (Human)	5018
Gene ID	5018
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 336-435 of human OXA1L (NP_005006.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15070
Protein Weight	49kDa
Shipping	Ice bag