

Anti-Oxgr1/Gpr99 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q96P68

Item Number: CM0002564



Introduction

Recombinant rabbit monoclonal antibody specific for human OXGR1/GPR99 (2-oxoglutarate receptor 1). Suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for studies on GPCR signaling, metabolism, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	OXGR1 / GPR99 Rabbit mAb
Remarks/Alias	aKGR; GPR80; GPR99; P2Y15; P2RY15; OXGR1 / GPR99
Species	Human
GeneID (Human)	27199
GeneID	27199
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human OXGR1 / GPR99 (Q96P68).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P68
Protein Weight	38 kDa
Shipping	Ice bag