


Anti-P2Rx2 Polyclonal Antibody For Wb, Elisa - Q9Ubl9

Item Number: CM0029409



Introduction

P2RX2 Rabbit pAb for WB and ELISA. Detects human, mouse, rat P2X purinoceptor 2, an ATP-gated ion channel involved in hearing and neurotransmission.

[Learn More](#)

Parameter	Value
Product Name	P2RX2 Rabbit pAb
Remarks/Alias	P2X2; DFNA41; P2RX2
Species	Human
Gene ID (Human)	22953
Gene ID	22953
Immunogen	Recombinant protein of human P2RX2 (amino acids 60-320)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBL9
Protein Weight	52kDa
Shipping	Ice bag