

Anti-P2Y12 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q9Cpv9

Item Number: CM0002872

Introduction

Anti-P2Y12 Rabbit Monoclonal Antibody validated for Western blot, IF-paraffin, and ELISA. Detects mouse and rat P2Y12, a platelet ADP receptor critical for hemostasis and thrombosis research. Ideal for coagulation studies.

[Learn More](#)

Parameter	Value
Product Name	P2Y12 Rabbit mAb
Remarks/Alias	P2Y12; 2900079B22Rik; 4921504D23Rik
Species	Mouse
GeneID (human)	28996
GeneID	70839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-347 of mouse P2Y12 (Q9CPV9)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9CPV9
Protein Weight	39kDa
Shipping	Ice bag