


Anti-P4Ha1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13674

Item Number: CM0013777



Introduction

Rabbit polyclonal antibody targeting human P4HA1. Validated for Western blot, IF/ICC, and ELISA. Detects a 61 kDa band; reacts with human, mouse, and rat. Essential for collagen biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	P4HA1 Rabbit pAb
Remarks/Alias	P4HA; P4HA1
Species	Human
GeneID (Human)	5033
GeneID	5033
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 235-534 of human P4HA1 (NP_001017962.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13674
Protein Weight	61kDa
Shipping	Ice bag