

Anti-P4Hb Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P07237

Item Number: CM0002132



Introduction

P4HB Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat P4HB, a 57 kDa protein disulfide-isomerase involved in redox regulation and protein folding.

[Learn More](#)

Parameter	Value
Product Name	P4HB Rabbit mAb
Remarks/Alias	DSI; GIT; PDI; PHDB; PDIA1; PO4DB; PO4HB; PROHB; CLCRP1; ERBA2L; P4Hbeta; P4HB
Species	Human
GeneID (Human)	5034
GeneID	5034
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 409-508 of human P4HB (P07237)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07237
Protein Weight	57kDa
Shipping	Ice bag