

Anti-P4Hb Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07237

Item Number: CM0016551



Introduction

Rabbit polyclonal antibody against human P4HB (Protein disulfide-isomerase) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studying protein folding, ER stress, and collagen biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	P4HB Rabbit pAb
Remarks/Alias	DSI; GIT; PDI; PHDB; PDIA1; PO4DB; PO4HB; PROHB; CLCRP1; ERBA2L; P4Hbeta; P4HB
Species	Human
GeneID	5034
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 329-508 of human P4HB (NP_000909.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07237
Protein Weight	57kDa
Shipping	Ice bag