

Anti-Padi2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2J8

Item Number: CM0015904



Introduction

PADI2 Rabbit pAb (Q9Y2J8) is a rabbit polyclonal antibody targeting protein-arginine deiminase type-2. Validated for WB and ELISA with human cross-reactivity, suitable for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	PADI2 Rabbit pAb
Remarks/Alias	PAD2; PDI2; PAD-H19; PADI2
Species	Human
Gene ID	11240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PADI2 (NP_031391.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y2J8
Protein Weight	76kDa
Shipping	Ice bag