

Anti-Padi4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Um07

Item Number: CM0022598



Introduction

Rabbit polyclonal antibody targeting human PADI4 (protein-arginine deiminase type-4). Validated in WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for studying chromatin decondensation, NETosis, and gene regulation.

[Learn More](#)

Parameter	Specification
Product Name	PADI4 Rabbit pAb
Remarks/Alias	PAD; PAD4; PDI4; PDI5; PADI5; PADI4
Species	Human
GeneID (Human)	23569
GeneID	23569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-663 of human PADI4 (NP_036519.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UM07
Protein Weight	74kDa
Shipping	Ice bag