


Anti-Pak3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O75914

Item Number: CM0009610



Introduction

PAK3 Rabbit mAb (O75914) validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for neuroscience research, targeting synaptic plasticity, dendritic spine morphology, and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PAK3 Rabbit mAb
Remarks/Alias	ARA; bPAK; MRX30; MRX47; OPHN3; PAK-3; XLID30; PAK3beta; beta-PAK; PAK3
Species	Human
GeneID (Human)	5063
GeneID	5063
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PAK3 (O75914)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75914
Protein Weight	62 kDa
Shipping	Ice bag