


Anti-Pan-Pkc Polyclonal Antibody For Wb, If/Icc, Elisa - P17252 / P05771 / Q05655 / Q02156 / P24723 / Q04759

Item Number: CM0022686



Introduction

High-quality rabbit polyclonal antibody against pan-PKC (PRKCA/PRKCB/PRKCD/PRKCE/PRKCH/PRKCQ), validated for WB, IF/ICC, ELISA in human, mouse, and rat. Ideal for PKC signaling research.

[Learn More](#)

Parameter	Value
Product Name	pan-PKC Rabbit pAb
Remarks/Alias	PRKCA/PRKCB/PRKCD/PRKCE/PRKCH/PRKCQ; pan-PKC
Species	Human
GeneID (human)	5578/5579/5580/5581/5583/5588
GeneID	5578/5579/5580/5581/5583/5588
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 572-672 of human PRKCA (NP_002728.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17252/P05771/Q05655/Q02156/P24723/Q04759
Protein Weight	59kDa/67kDa/74- 83kDa
Shipping	Ice bag