

Anti-Pank2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bz23

Item Number: CM0021834



Introduction

Rabbit polyclonal antibody against human PANK2, for WB and ELISA. Cross-reacts with mouse. Immunogen: synthetic peptide aa 300-400. Useful for coenzyme A synthesis and neurodegeneration research.

[Learn More](#)

Parameter	Specification
Product Name	PANK2 Rabbit pAb
Remarks/Alias	HSS; HARP; PKAN; NBIA1; C20orf48; PANK2
Species	Human
GeneID (Human)	80025
GeneID	80025
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human PANK2 (NP_705902.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9BZ23
Protein Weight	63kDa
Shipping	Ice bag