

Anti-Parl Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H300

Item Number: CM0030841

Introduction

Anti-PARL rabbit polyclonal antibody for WB, IF/ICC, ELISA. Recognizes human PARL, a mitochondrial protease involved in apoptosis regulation. Ideal for apoptosis and mitochondrial function studies.

[Learn More](#)

Parameter	Value
Product Name	PARL Rabbit pAb
Remarks/Alias	PSARL; PSARL1; RHBDL1; PRO2207; PSENIP2; PARL
Species	Human
GeneID (Human)	55486
GeneID	55486
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human PARL (NP_001032728.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9H300
Protein Weight	42kDa
Shipping	Ice bag