

Anti-Pax5 Polyclonal Antibody For Wb, Elisa - Q02548

Item Number: CM0017032



Introduction

High-quality rabbit polyclonal antibody against human PAX5 (BSAP), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for B-cell development and EBV research.

[Learn More](#)

Parameter	Value
Product Name	PAX5 Rabbit pAb
Remarks/Alias	ALL3; BSAP; PAX-5; PAX5
Species	Human
GeneID (Human)	5079
GeneID	5079
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PAX5 (NP_057953.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02548
Protein Weight	42kDa
Shipping	Ice bag