

Anti-Pax5 Polyclonal Antibody For Wb, Ip, Elisa - Q02548

Item Number: CM0026478



Introduction

Rabbit polyclonal antibody to human PAX5 transcription factor, suitable for WB, IP, ELISA. Detects a 42kDa band. Target essential for B-cell development and EBV genome maintenance.

[Learn More](#)

Parameter	Value
Product Name	PAX5 Rabbit pAb
Remarks/Alias	ALL3; BSAP; PAX-5; PAX5
Species	Human
GeneID (Human)	5079
GeneID	5079
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-300 of human PAX5 (NP_057953.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q02548
Protein Weight	42kDa
Shipping	Ice bag