

Anti-Pcdha12 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Un75

Item Number: CM0030073

Introduction

This rabbit polyclonal antibody targets protocadherin alpha-12 (PCDHA12), a neuronal cell adhesion protein. Validated for WB and ELISA, it recognizes human, mouse, and rat PCDHA12 (~102kDa). Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PCDHA12 Rabbit pAb
Remarks/Alias	PCDH-ALPHA12; PCDHA12
Species	Human
GeneID (Human)	56137
GeneID	56137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 713-792 of human PCDHA12 (NP_114070.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UN75
Protein Weight	102kDa
Shipping	Ice bag