

Anti-Pcsk1 Polyclonal Antibody For Wb, If/Icc, Elisa - P29120

Item Number: CM0018529



Introduction

Affinity-purified rabbit polyclonal antibody against human PCSK1 (Neuroendocrine convertase 1). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying prohormone processing and neuroendocrine disorders.

[Learn More](#)

Parameter	Value
Product Name	PCSK1 Rabbit pAb
Remarks/Alias	PC1; PC3; NEC1; SPC3; PC1/3; BMIQ12; PCSK1
Species	Human
GeneID (human)	5122
GeneID	5122
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-170 of human PCSK1 (NP_000430.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P29120
Protein Weight	84kDa
Shipping	Ice bag