

Anti-PCYT1A Monoclonal Antibody For Wb, Elisa - P49585

Item Number: CM0004469

Introduction

Anti-PCYT1A Rabbit Monoclonal Antibody for WB and ELISA, validated in human samples. This antibody targets choline-phosphate cytidylyltransferase A, a key enzyme in phosphatidylcholine synthesis, suitable for metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	PCYT1A Rabbit mAb
Remarks/Alias	CT; CTA; CCTA; CTPCT; PCYT1; SMDCRD; CCTalpha; PCYT1A
Species	Human
GeneID (Human)	5130
GeneID	5130
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 268-367 of human PCYT1A (P49585)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49585
Protein Weight	42kDa
Shipping	Ice bag