

Anti-Pd-L1/Cd274 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Nzq7

Item Number: CM0022915



Introduction

High-quality anti-PD-L1/CD274 rabbit polyclonal antibody validated for WB, IF-P, IHC-P, ELISA. Detects human, mouse, rat PD-L1; suitable for cancer immunotherapy and immune tolerance studies.

[Learn More](#)

Parameter	Value
Product Name	PD-L1/CD274 Rabbit pAb
Remarks/Alias	B7-H; B7H1; PDL1; PD-L1; hPD-L1; PDCD1L1; PDCD1LG1; PD-L1/CD274
Species	Human
GeneID (Human)	29126
GeneID	29126
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 190-290 of human PD-L1/CD274 (NP_054862.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZQ7
Protein Weight	33kDa
Shipping	Ice bag