


Anti-Pde10A Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y233

Item Number: CM0004398



Introduction

Rabbit monoclonal antibody targeting human PDE10A, validated for WB and ELISA. Cross-reacts with mouse and rat. PDE10A regulates cAMP/cGMP signaling, highly expressed in striatum. Ideal for neuroscience and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	PDE10A Rabbit mAb
Remarks/Alias	ADSD2; IOLOD; HSPDE10A; PDE10A19; LINC00473; PDE10A
Species	Human
Gene ID	10846
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 439-779 of human PDE10A (NP_006652.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y233
Protein Weight	88kDa
Shipping	Ice bag