


Anti-Pde3B Rabbit Polyclonal Antibody For Wb, Elisa - Q13370

Item Number: CM0025800



Introduction

Affinity-purified rabbit anti-PDE3B polyclonal antibody validated for WB and ELISA. Detects human and mouse PDE3B, a dual-specificity phosphodiesterase involved in angiogenesis and cardiac function. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	PDE3B Rabbit pAb
Remarks/Alias	HcGIP1; cGIPDE1; PDE3B
Species	Human
Gene ID (Human)	5140
Gene ID	5140
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-660 of human PDE3B (NP_000913.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13370
Protein Weight	124kDa
Shipping	Ice bag