

Anti-Pdgfra Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P16234

Item Number: CM0010439



Introduction

Rabbit monoclonal PDGFRa antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, rat PDGFRa (CD140A). Recombinant immunogen (aa 24-528). Ideal for cancer and developmental research.

[Learn More](#)

Parameter	Value
Product Name	PDGFRa Rabbit mAb
Remarks/Alias	CD140A, PDGFR2, PDGFR-2
Species	Human
GeneID (Human)	5156
GeneID	5156
Immunogen	Recombinant protein (24-528aa of human PDGFRa)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16234
Protein Weight	123kDa
Shipping	Ice bag