


Anti-Pdlim1/Clp36 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00151

Item Number: CM0025434



Introduction

Rabbit polyclonal antibody against human PDLIM1/CLP36, validated for WB, IF/ICC, and ELISA. Detects a 36kDa cytoskeletal adapter protein involved in stress fiber organization. Cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	PDLIM1/CLP36 Rabbit pAb
Remarks/Alias	CLIM1; CLP36; CLP-36; hCLIM1; HEL-S-112; PDLIM1/CLP36
Species	Human
Gene ID (Human)	9124
Gene ID	9124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-329 of human PDLIM1/CLP36 (NP_066272.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00151
Protein Weight	36kDa
Shipping	Ice bag