

Anti-Pdxdc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6P996

Item Number: CM0019876

Introduction

High-quality rabbit polyclonal antibody targeting human PDXDC1, a pyridoxal-dependent decarboxylase. Validated for WB, IHC-P, and ELISA; cross-reacts with human, mouse, and rat. 87 kDa.

[Learn More](#)

Parameter	Value
Product Name	PDXDC1 Rabbit pAb
Remarks/Alias	LP8165; PDXDC1
Species	Human
GeneID (Human)	23042
GeneID	23042
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 489-788 of human PDXDC1 (NP_055842.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P996
Protein Weight	87kDa
Shipping	Ice bag