

Anti-Pepck/Pck2 Polyclonal Antibody For Wb And Elisa - Q16822

Item Number: CM0025666

Introduction

Rabbit polyclonal antibody against mitochondrial PEPCK (PCK2), validated for WB and ELISA in human, mouse, and rat. Ideal for gluconeogenesis and ferroptosis research.

[Learn More](#)

Parameter	Value
Product Name	PEPCK/PCK2 Rabbit pAb
Remarks/Alias	PEPCK; PEPCK2; PEPCK-M; PEPCK/PCK2
Species	Human
GeneID	5106
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PEPCK/PCK2 (NP_004554.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16822
Protein Weight	71kDa
Shipping	Ice bag