

Anti-Pepck/Pck2 Rabbit Monoclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - Q16822

Item Number: CM0004921

Introduction

Recombinant anti-PEPCK/PCK2 rabbit monoclonal antibody, KO validated, for western blot, immunofluorescence, and ELISA. Detects human, mouse, and rat mitochondrial PCK2. Ideal for gluconeogenesis and ferroptosis research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PEPCK/PCK2 Rabbit mAb
Remarks/Alias	PEPCK; PEPCK2; PEPCK-M; K2
Species	Human
GeneID (Human)	5106
GeneID	5106
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-640 of human PEPCK/PCK2 (Q16822)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16822
Protein Weight	71kDa
Shipping	Ice bag