

Anti-Pepd Rabbit Polyclonal Antibody For Wb, Elisa - P12955

Item Number: CM0028311



Introduction

High-quality rabbit polyclonal antibody against human PEPD (prolidase). Validated for Western blotting and ELISA applications. Cross-reacts with mouse and rat. Ideal for investigating collagen metabolism and related research. Approximate molecular weight: 54kDa.

[Learn More](#)

Parameter	Value
Product Name	PEPD Rabbit pAb
Remarks/Alias	PROLIDASE; PEPD
Species	Human
Gene ID (Human)	5184
Gene ID	5184
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 288-424 of human PEPD (NP_000276.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12955
Protein Weight	54kDa
Shipping	Ice bag