

Anti-Pepsinogen Ii (Pg Ii) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P20142

Item Number: CM0026719



Introduction

Anti-Pepsinogen II (PG II) rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat gastricsin. Ideal for pepsinogen C studies in digestion and protease research.

[Learn More](#)

Parameter	Value
Product Name	Pepsinogen II (PG II) Rabbit pAb
Remarks/Alias	PEPC; PGII; Pepsinogen II (PG II)
Species	Human
GeneID (human)	5225
GeneID	5225
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 73-172 of human Pepsinogen II (PG II) (NP_001159896.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20142
Protein Weight	42kDa
Shipping	Ice bag