

Anti-Per2 Polyclonal Antibody For Wb, Ihc-P, Elisa, Chip - O15055

Item Number: CM0016411



Introduction

Rabbit polyclonal antibody against human PER2 (Period circadian protein homolog 2), validated for WB, IHC-P, ELISA, and ChIP. Cross-reacts with human, mouse, and rat. Suitable for circadian rhythm research.

[Learn More](#)

Parameter	Value
Product Name	PER2 Rabbit pAb
Remarks/Alias	FASPS; FASPS1; PER2
Species	Human
Gene ID (human)	8864
Gene ID	8864
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 390-565 of human PER2 (NP_073728.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	O15055
Protein Weight	137kDa
Shipping	Ice bag