


Anti-Peroxiredoxin 4 (Prdx4) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13162

Item Number: CM0022599



Introduction

Rabbit polyclonal antibody against human Peroxiredoxin 4 (PRDX4). Validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat PRDX4. Suitable for oxidative stress and redox signaling research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 4 (PRDX4) Rabbit pAb
Remarks/Alias	PRX-4; AOE372; AOE37-2; HEL-S-97n; Peroxiredoxin 4 (PRDX4)
Species	Human
GeneID (human)	10549
GeneID	10549
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 10-110 of human Peroxiredoxin 4 (PRDX4) (NP_006397.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13162
Protein Weight	31kDa
Shipping	Ice bag