


Anti-Peroxiredoxin 5 (Prdx5) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P30044

Item Number: CM0014246



Introduction

High-quality Anti-Peroxiredoxin 5 (PRDX5) Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects Human, Mouse, Rat PRDX5. Ideal for oxidative stress, antioxidant defense, and peroxisomal research.

[Learn More](#)

Parameter	Value
Product Name	Peroxiredoxin 5 (PRDX5) Rabbit pAb
Remarks/Alias	PLP; ACR1; B166; PRXV; PMP20; PRDX6; prx-V; SBB110; AOEB166; HEL-S-55; Peroxiredoxin 5 (PRDX5)
Species	Human
GeneID (Human)	25824
GeneID	25824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 53-214 of human Peroxiredoxin 5 (PRDX5) (NP_036226.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30044
Protein Weight	22kDa
Shipping	Ice bag