

Anti-Pex10 Rabbit Polyclonal Antibody For Wb And Elisa - O60683

Item Number: CM0021139

Introduction

This rabbit polyclonal antibody targets human PEX10, validated for WB and ELISA with cross-reactivity to mouse. Suitable for investigating peroxisome biogenesis and PEX5 recycling pathways in research.

[Learn More](#)

Parameter	Value
Product Name	PEX10 Rabbit pAb
Remarks/Alias	NALD; PBD6A; PBD6B; RNF69; PEX10
Species	Human
Gene ID (Human)	5192
Gene ID	5192
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 10-160 of human PEX10 (NP_002608.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60683
Protein Weight	37kDa
Shipping	Ice bag