


Anti-Pex19 Rabbit Monoclonal Antibody For Wb And Elisa - P40855

Item Number: CM0002161



Introduction

CamelBio's anti-PEX19 rabbit monoclonal antibody, validated for western blotting and ELISA, specifically recognizes human PEX19. Essential for peroxisomal biogenesis studies, cell cycle research, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PEX19 Rabbit mAb
Remarks/Alias	PXF; HK33; PMP1; PMPI; PXMP1; PBD12A; D152223E; PEX19
Species	Human
Gene ID	5824
Immunogen	Synthetic peptide within amino acids 1-100 of human PEX19 (P40855)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P40855
Protein Weight	33kDa
Shipping	Ice bag