

Anti-Pgbd1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Js3

Item Number: CM0023551

Introduction

Rabbit polyclonal antibody against human PGBD1 (SCAND4), validated for WB and ELISA. PGBD1 is a nuclear transposase-derived protein that regulates neuronal gene expression and paraspeckle assembly. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	PGBD1 Rabbit pAb
Remarks/Alias	SCAND4; HUCEP-4; dJ874C20.4; PGBD1
Species	Human
GeneID (human)	84547
GeneID	84547
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 150-450 of human PGBD1 (NP_115896.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96JS3
Protein Weight	93kDa
Shipping	Ice bag